

1  
2  
3  
4

**FINAL REPORT**  
**PDX-LIB Labor Saver**  
**AOAC-RI Performance Tested Methods**

5 **Abstract**

6 PDX-LIB Labor Saver method is derived from the previously approved PDX-LIB  
7 method (AOAC 040501) based on the fact that the industry wanted a “simpler” method.  
8 Derivation provided elimination of laborious handling and increased sampling size.  
9 Paradigm Diagnostics’ *Listeria* Indicator Broth (PDX-LIB) is a specialty medium  
10 intended to be used as a screening tool for the detection of *Listeria* bacteria from  
11 environmental samples collected from food processing facilities. PDX-LIB is distinct  
12 because it is a single stage enrichment that can be used as a pre-screener for confirmatory  
13 tests. Interpretation of results is based on visible color change (darkening from yellow  
14 amber to brown black in color) after 30 to 48 hours of incubation at 37 °C. Performance  
15 of PDX-LIB Labor Saver has been tested in five different categories; inclusivity-  
16 exclusivity, method comparison, ruggedness lot to lot variability and shelf life. There  
17 were 67 *Listeria* and 37 non-*Listeria* species tested in the inclusivity-exclusivity studies.  
18 Inclusivity-exclusivity results suggested that all 67 species of *Listeria* tested gave  
19 positive results at less than 52 CFU. Only a few *Enterococcus* species gave false positive  
20 results at high cell counts. Method comparison studies were done for *Listeria*  
21 *monocytogenes* (Lm) in the presence of 1 log excess *E.coli* on sealed concrete, *L. innocua*  
22 (Li) on ceramic tile *L. ivanovii* (Liv) on stainless steel, and *L. welshmeri* (Lw) on plastic.  
23 Results of method comparison studies suggested that PDX-LIB Labor Saver performed  
24 better than the USDA method. Results of the ruggedness studies suggested deviations  
25 from the test parameters did not affect the test performance significantly. Lot-to-lot  
26 variability studies showed that there were no performance differences among three  
27 different production lots. PDX-LIB has been shown to pass QC criteria after 9 months of  
28 storage at refrigerated conditions.

29  
30 **Method Authors**

31 H. Cem Yurttas\*, Jessica L. Saylor, Leena Griffith, Karen Larson, Alan Olstein and  
32 Joellen Feirtag

33 \*Corresponding author  
34

**Submitting Company**

Research and Development Team  
Paradigm Diagnostics LLC.  
1334 Eckles Ave  
St. Paul, MN 55108

35

1

## 2 **1 Scope of Method**

3

4 PDX-LIB Labor Saver is a presumptive *Listeria* test that has been specifically designed  
5 to monitor presence of *Listeria* in environmental samples from the food processing  
6 environment. The test contains a liquid media formulated with antibiotics, growth  
7 enhancers, and color changing compounds specific to the growth of *Listeria* species.  
8 PDX-LIB Labor Saver has been vigorously tested not only at low levels of pure *Listeria*  
9 cultures (1-100 CFU/mL) but also in naturally contaminated environmental samples from  
10 the food industry during the initial developmental studies. Applicability of PDX-LIB  
11 Labor Saver is limited for selected common *Listeria* spp (Lm, Li, Liv, and Lw) on  
12 selected common surface types (sealed concrete, ceramic tile, stainless steel and plastic).  
13 In order to assess the performance characteristics, PDX-LIB has been tested in five  
14 different performance categories; 1) Inclusivity and exclusivity studies 2) Method  
15 comparison studies 3) Ruggedness studies 4) Lot to lot variability and 5) Shelf stability  
16 studies.

17

## 18 **2 Definitions**

19

20 2.1 Chi-squared ( $X^2$ ) Value: The Chi-Square value for unpaired samples is frequently  
21 used in Official Methods to determine if two methods are equivalent. The formula  
22 is;

23

$$24 X^2 = N * \{[(N_{12})(N_{21}) - (N_{11})(N_{22})]^2 / N\} / \{(N_{12} + N_{11})(N_{12} + N_{22})(N_{22} + N_{21})(N_{11} + N_{21})\}$$

25

26  $N_{12}$  = samples positive by Test method

27  $N_{11}$  = samples positive by the reference method

28  $N_{22}$  = samples negative by Test method

29  $N_{21}$  = samples negative by the reference method

30

## 31 **3 Principle:**

32

33 The principle of the PDX-LIB is based on colorometric changes in enrichment medium  
34 in the presence of *Listeria* spp utilizing a unique blend of antibiotics, growth enhancers  
35 and *Listeria* specific color indicators working synergistically. A light brown to black  
36 color formation within 30 hours indicates presence of *Listeria* spp in an environmental  
37 samples. PDX-LIB has a limit of detection for most *Listeria* spp between (>1-100  
38 CFU/mL) within 30 hours.

39

### 40 **3.1 General Information**

41 *Listeria* species are naturally found in the environment including animals, plants, soil,  
42 dust, and silage. *Listeria* spp. are gram-positive non-spore forming rods and have tumble  
43 motility. There are six different species of *Listeria*; *Listeria monocytogenes*, (Lm) *L.*  
44 *innocua* (Li), *L. ivanovii* (Liv), *L. welshmeri* (Lw), *L. grayi* (Lg), and *L. seeligeri* (Ls). Of  
45 these six species, *L. monocytogenes* and *L. ivanovii* are considered pathogens. *L.*  
46 *monocytogenes* is the most deadly of the pathogenic species to humans. *L. ivanovii* was

1 once considered to be a separate serotype of *L. monocytogenes* until it was reclassified in  
2 1984 as a new species (1). Because of this close connection to the human pathogen, *L.*  
3 *ivanovii* although more likely to cause sickness in animals, can be considered a rare but  
4 viable threat to humans (2).

5  
6 For a non-spore forming bacterium, *L. monocytogenes* is relatively resistant to heating  
7 and drying and can multiply at temperatures as low as 3°C (3). Although it is able to  
8 withstand some heating and drying, it can still be injured during the course of sublethal  
9 heating, freezing, and exposure to detergents. It is the detection of this type of *Listeria*  
10 that is particularly important. Sub-lethally injured organisms can undergo repair in food  
11 product and cause illness to consumers (4).

12  
13 *Listeria spp.* can typically be found in uncooked meats and vegetables, un-pasteurized  
14 milk, ice cream, soft cheeses like Swiss and Mexican, and ready to eat food products  
15 (RTE) (5). Animals carry these bacteria in their intestinal tract without showing signs of  
16 illness. These bacteria can then be transferred to meat as well as dairy products during  
17 food processing.

18  
19 *Listeriosis* is the disease caused by consuming food contaminated with *L.*  
20 *monocytogenes*. It is one of the leading causes of death not because it occurs often but  
21 because of the extremely high mortality rate (6). In the United States alone, there are  
22 over 2,500 reported cases of *Listeriosis* every year; of these there is approximately 500  
23 deaths. Individuals who are immuno-compromised, young, old, or pregnant should be  
24 careful not to consume adulterated products. Symptoms may include fever, muscle  
25 aches, and diarrhea. If the nervous system becomes infected symptoms may include  
26 headache, stiff neck, confusion, loss of balance and convulsions. In pregnant women,  
27 *listeriosis* can produce flu like illness and can lead to miscarriage or stillbirth of the fetus  
28 (5). Because of the seriousness of the illness caused by *Listeria* contamination,  
29 government agencies have taken steps to reduce incidence of *Listeria* outbreaks including  
30 a “zero-tolerance” policy for *Listeria* in 25g samples of food product (1). However,  
31 several European countries instituted an action level for 100 CFU/g for *L.monocytogenes*  
32 (7). Furthermore according to Chen *et al*, 2003(8), foods containing low levels (<100  
33 cfu/g) of *L.monocytogenes* pose very little risk.

34  
35 Between January 2004 and April of 2004, there has been six recalls on various food  
36 products contaminated with *Listeria* (9). Each of these recalls cost companies a  
37 tremendous amount of money as well as reputation among consumers.

38  
39 Current methods for enrichment and detection of *Listeria spp.* can take 3-5 days for  
40 presumptive results and an additional 2-3 days for confirmation. This method is time  
41 consuming and requires a great amount of work (3). It is also important to note that the  
42 detection of *Listeria* is limited by the abilities of the method. Detection and enrichment  
43 medias that rely heavily on selective agents may not be able to culture out sub-lethally  
44 injured *Listeria* organisms. Although these organisms occasionally do not show up in the  
45 culture method, they can persist in the food product and cause an outbreak to occur (4).

1 It is clear that a more rapid method with equal or better sensitivity is needed to control  
2 *Listeria* in the food industry.

3

#### 4 **4 Test Kit Information**

5 4.1 Kit Name: PDX-LIB

6 4.2 Catalog Numbers: 25003-25, 25004-100, 25005-100, 25009-50

7 4.2.1 In the US and Canada

##### 8 **Paradigm Diagnostics, LLC**

9 1360 University Ave. W, Suite 455

10 St. Paul, MN 55104

11 Web: [www.pdx-inc.com](http://www.pdx-inc.com)

12

##### 13 **Hardy Diagnostics**

14 Christopher A. Catani MT(ASCP), RM(ACM)

15 Director of Sales, Marketing & Technical Support

16

17 Ph: (800) 266-2222 Ext. 5696

18 Fx: (805) 614-9274

19 [CataniC@HardyDiagnostics.com](mailto:CataniC@HardyDiagnostics.com)

20 [www.HardyDiagnostics.com](http://www.HardyDiagnostics.com)

21

22 Hardy Diagnostics

23 1430 West McCoy Lane

24 Santa Maria, CA 93455

25 4.2.2 New Zealand

##### 26 **Food Tech Solutions Ltd**

27

28 Peter J. Bailey

29 Director of Sales and Marketing,

30 Technical Support

31

32 Ph: (64) 9 576 7326

33 Fax: (64) 9 576 7005

34 [foodtech@ihug.co.nz](mailto:foodtech@ihug.co.nz)

35 [www.fts.co.nz](http://www.fts.co.nz)

36

#### 37 **4.3 Test Kit Reagents**

38 4.3.1 Aseptically filled 20 mL of PDX-LIB tubes

39

#### 40 **5 Additional Supplies and Reagents**

41 5.1 Tecra Enviroswabs

42 5.2 MOX plates: Modified Oxford Agar plates for optional confirmation step

43

## 6 Apparatus

6.1 Incubator at 37 °C (not included in the test kit)

## 7 Standard Reference Materials

All reference materials required by USDA document MLG 8.04 have been used as necessary (Source: [http://www.fsis.usda.gov/Ophs/Microlab/Mlg\\_8\\_04.pdf](http://www.fsis.usda.gov/Ophs/Microlab/Mlg_8_04.pdf))

## 8 Standard Solutions

Not applicable

## 9 Safety Precautions

PDX-LIB Labor Saver test should be done by personnel with basic microbiology laboratory training. Material Safety Data Sheets are available upon request. All materials used should be handled and disposed of as potentially infectious material. Autoclaving is the preferred method of disposal. If autoclaving is not available, disinfectant solutions should be used. *Listeria monocytogenes* is a dangerous human pathogen. When handling samples that possibly contain *L.monocytogenes*, extreme care should be taken to contain the samples and the enriched samples (presumptive positive tubes). Immuno-compromised individuals and pregnant women are particularly endangered by exposure to *L.monocytogenes* and should not be allowed in the vicinity of the testing.

## 10 General Preparation

The pure cultures were obtained from R-Tech Labs, Arden Hills, MN, American Typed Culture Collection (ATCC), Manassas, VA, Veterinary Medical Diagnostics Labs at University of Minnesota, St. Paul, MN, and Microbiology Labs of Dr. Diez and Dr. Feirtag at the University of Minnesota, St. Paul, MN. Original pure cultures were grown overnight in sterile Tryptic Soy Broth (TSB) at 32-37 °C. In order to have reproducible cultures, 50 % glycerol stock solutions of cultures were prepared and stored in a freezer. Fifty percent glycerol stock cultures were prepared by diluting 500 µL of pure culture grown overnight with sterile 500 µL glycerol. Then 200 uL portions were filled into sterile 2 mL centrifuge tubes, capped and kept in a freezer at -15 °C until the day before use. On the day before use, a loop-full of freezer-stored stock cultures were inoculated into 5 mL of sterile TSB, the same loop was streaked onto tryptic soy agar (TSA) containing ferric ammonium citrate and esculin (TSAIE). The TSB and TSAIE were incubated overnight at 37 °C. TSAIE plates were checked for contamination, based on β-glucosidase activity and colony morphology. If the TSAIE plates suggested there was no apparent contamination, logarithmic level dilutions were made into sterile peptone solutions (0.1 % peptone in distilled water). Aliquots containing 100 µL of 6 and 7 log dilutions from overnight cultures were plated onto TSAIE for estimating the cell concentration (CFU/mL) and checked for contamination.

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46

## **11 Sample Preparation**

### **11.1 Sampling for AOAC-RI Guideline Studies:**

Dilutions of pure cultures aforementioned were used for all AOAC-RI recommended studies, which were inclusivity-exclusivity, method comparison, ruggedness and lot to-lot variability and shelf stability studies.

## **12 Analysis**

- Take an environmental sample following USDA recommendations for sampling. Return the sampling device back into the original sterile container.
- Aseptically add one unit (20 mL) of PDX-LIB media to the sterile container to fully submerge the applicator tip.
- Incubate in upright position for 30 to 48 hours at 37 °C.
- Check the color of the media. Sample is called presumptive positive for *Listeria* species if color of the media changed from yellow/amber color to black/brown color.

## **13 Interpretation and Test Result Report**

Light to dark brown-black colored media at the end of 30 to 48 hour indicates presumptive presence of *Listeria spp.* In order to confirm the negatives, all yellow colored (presumptive negative) tubes have to be incubated for total of 48 hours. It is recommended that at least one negative control be run in each set of analysis. A negative control is an unused sampling device containing one unit of PDX-LIB, incubated alongside the samples.

## **14 Internal Validation Studies**

### **14.1. Inclusivity –Exclusivity Studies**

Serial dilutions from overnight grown pure cultures of 50 different *Listeria* and 30 different non-*Listeria* cultures were made into sterile peptone. One hundred microliters from log –6 dilution (~1.0E+3 CFU/mL) were plated onto TSAIE plates for estimating the number of cells tested. For the inclusivity studies 0.5 mL (estimated cell concentration ranging from 1 to 52 CFU of *Listeria*) of a dilution aliquot was aseptically transferred onto the tip of the sampling sponge. For the exclusivity studies 0.5 mL (estimated cell concentration ranging from 10,000 to 100,000 CFU of non-*Listeria*) of a dilution aliquot was aseptically transferred onto the tip of the sampling sponge. Then sampling device was placed back into the sterile tube and submerged in 20 mL of PDX-LIB media. The results are given in Table 1a and 1 b respectively.

1 Table 1 a. Results of Inclusivity Tests for PDX-LIB

Description and code	Source	Color	Presumptive Result	API
<i>Listeria innocua</i> 33090	Cow brain	Black	+	Li
<i>Listeria welshmeri</i> 43551 (B)	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 19112 (B)	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 3528	Food/Food Environment	Black	+	Lm
<i>Listeria innocua</i> 2242	Food/Food Environment	Black	+	Li
<i>Listeria monocytogenes</i> 3522	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2421	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2397	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 4b rt 652 (B)	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2395	Food/Food Environment	Black	+	Lm
<i>Listeria ivanovii</i> 700402	Food/Food Environment	Black	+	Liv
<i>Listeria monocytogenes</i> 2392	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2415	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 19118 (B)	Food/Food Environment	Black	+	Lm
<i>Listeria grayi</i> rt721	Food/Food Environment	Black	+	Lg
<i>Listeria monocytogenes</i> 4b 19115	Human	Black	+	Lm
<i>Listeria grayi</i> 700545 Lab	Unspecified	Black	+	Lg
<i>Listeria monocytogenes</i> 1/2 a rt 651 (B)	Food/Food Environment	Black	+	Lm
PDX-Lm47	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2417	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 1/2 b rt 541 (B)	Food/Food Environment	Black	+	Lm
PDX-Lm48	Food/Food Environment	Black	+	Lm
PDX-Lm44	Food/Food Environment	Black	+	Lm
PDX-Ls4	Unknown	Black	+	Ls
<i>Listeria welshmeri</i> 35967 Lab	Unknown	Black	+	Lw
<i>Listeria monocytogenes</i> 2422	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2413	Food/Food Environment	Black	+	Lm
U of MN-Vet Med	Unspecified	Black	+	Lm
U of MN-Vet Med	Unspecified	Black	+	Lm
<i>Listeria innocua</i> 2241	Food/Food Environment	Black	+	Li
<i>Listeria monocytogenes</i> 2426	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2578	Food/Food Environment	Black	+	Lm
PDX-Li10	Food/Food Environment	Black	+	Li
<i>Listeria monocytogenes</i> 2388	Food/Food Environment	Black	+	Lm
<i>Listeria seeligeri</i> 2232	Unspecified	Black	+	Ls
<i>Listeria monocytogenes</i> 2427	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 3742	Food/Food Environment	Black	+	Lm
U of MN-Vet Med	Unspecified	Black	+	Lm
<i>Listeria innocua</i> 51742 Lab	Unknown	Black	+	Li
PDX-Lw3	Unknown	Black	+	Lw
U of MN-Vet Med	Unspecified	Black	+	Lm
U of MN-Vet Med	Unspecified	Black	+	Lm
PDX-Lm46	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 19114	Food	Black	+	Lm
<i>Listeria monocytogenes</i> 19117	Food/Food Environment	Black	+	Lm
PDX-Lm43	Food/Food Environment	Black	+	Lm
PDX-Lm45	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2389	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2396	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 2424	Food/Food Environment	Black	+	Lm
<i>Listeria monocytogenes</i> 15313	Food/Food Environment	Black	+	Lm
U of MN-Vet Med	Unspecified	Black	+	Lm
<i>Listeria monocytogenes</i> 1914	Food/Food Environment	Black	+	Lm
<i>Listeria grayi</i> 19120	Chincilla feces	Black	+	Lg
<i>Listeria monocytogenes</i> 2349	Food/Food Environment	Black	+	Lm
U of MN-Vet Med	Unspecified	Black	+	Lm
<i>Listeria innocua</i> 2249	Food/Food Environment	Black	+	Li
<i>Listeria monocytogenes</i> 3550	Food/Food Environment	Black	+	Lm
U of MN-Vet Med	Unspecified	Black	+	Lm
<i>Listeria monocytogenes</i> 2410	Food/Food Environment	Black	+	Lm
PDX-Li 2	Unknown	Black	+	Li
<i>Listeria monocytogenes</i> 2404	Food/Food Environment	Black	+	Lm
<i>Listeria innocua</i> 3757	Food/Food Environment	Black	+	Li
U of MN-Vet Med	Unspecified	Black	+	Lm
<i>Listeria innocua</i> 3181	Food/Food Environment	Black	+	Li
<i>Listeria innocua</i> 3254	Food/Food Environment	Black	+	Li
<i>Listeria ivanovii</i> 49954	Food, France	Black	+	Liv

2

1

2 Table 1 b Results of Exclusivity Tests for PDX-LIB

Description and code	Source	Color at 30 hours	Color at 48 hours	Presumptive Result	Confirmation
<i>E.coli</i> O157:H7 43895	ATCC-undefined	Yellow	Yellow	-	Non-Listeria
<i>E. coli</i> 25922	Clinical	Yellow	Yellow	-	Non-Listeria
<i>E. coli</i> 10798	Feces	Yellow	Yellow	-	Non-Listeria
<i>E. coli</i> 51739	ATCC-undefined	Yellow	Yellow	-	Non-Listeria
<i>Bacillus cereus</i> 11778	ATCC-undefined	Yellow	Yellow	-	Non-Listeria
<i>Bacillus cereus</i>	Feirtag Culture Collection	Yellow	Yellow	-	Non-Listeria
<i>Bacillus cereus</i>	Feirtag Culture Collection	Yellow	Yellow	-	Non-Listeria
<i>Bacillus stearothermophilus</i>	Feirtag Culture Collection	Yellow	Yellow	-	Non-Listeria
<i>Bacillus stearothermophilus</i>	Feirtag Culture Collection	Yellow	Yellow	-	Non-Listeria
<i>Bacillus subtilis</i>	Feirtag Culture Collection	Yellow	Yellow	-	Non-Listeria
<i>Bacillus circulans</i>	Feirtag Culture Collection	Yellow	Yellow	-	Non-Listeria
<i>Bacillus licheniformis</i> 12759	Cassava tuber	Yellow	Yellow	-	Non-Listeria
<i>Bacillus licheniformis</i> 14580	ATCC-undefined	Yellow	Yellow	-	Non-Listeria
<b><i>Enterococcus hire</i> 8043</b>	<b>ATCC-undefined</b>	<b>Black</b>	<b>Black</b>	<b>+</b>	Non-Listeria
<b><i>Enterococcus hire</i> 49611</b>	<b>Pig out</b>	<b>Black</b>	<b>Black</b>	<b>+</b>	Non-Listeria
<b><i>Enterococcus hire</i> 49612</b>	<b>ATCC-undefined</b>	<b>Black</b>	<b>Black</b>	<b>+</b>	Non-Listeria
<b><i>Enterococcus hire</i> 35220</b>	<b>Cattle dung</b>	<b>Black</b>	<b>Black</b>	<b>+</b>	Non-Listeria
<i>Enterobacter aerogenes</i> 13048	Sputum	Yellow	Yellow	-	Non-Listeria
<i>Enterobacter cloacae</i> 13047	Spinal fluid	Yellow	Yellow	-	Non-Listeria
<b><i>Enterococcus faecalis</i> 29212</b>	<b>Urine</b>	<b>Yellow</b>	<b>Black</b>	<b>+</b>	Non-Listeria
<i>Enterococcus faecalis</i> 19433	ATCC-undefined	Yellow	Yellow	-	Non-Listeria
<i>Enterococcus durans</i> 19432	ATCC-undefined	Yellow	Yellow	-	Non-Listeria
<b><i>Enterococcus gallinarum</i> 49573</b>	<b>Chicken intestine</b>	<b>Yellow</b>	<b>Black</b>	<b>+</b>	Non-Listeria
<b><i>Enterococcus avium</i> 14025</b>	<b>ATCC-undefined</b>	<b>Black</b>	<b>Black</b>	<b>+</b>	Non-Listeria
<i>Enterococcus casseliflavus</i> 25788	Plants	Yellow	Yellow	-	Non-Listeria
<i>Staphylococcus aureus</i> 36548	ATCC-undefined	Yellow	Yellow	-	Non-Listeria
<i>Staphylococcus aureus</i> 8095	R-tech labs-undefined	Yellow	Yellow	-	Non-Listeria
<i>Staphylococcus aureus</i> 25923	Clinical	Yellow	Yellow	-	Non-Listeria
<i>Staphylococcus aureus</i> 51740	Margarine	Yellow	Yellow	-	Non-Listeria
<i>Klebsiella pneumoniae</i> 13883	ATCC-undefined	Yellow	Yellow	-	Non-Listeria
<i>Klebsiella pneumoniae</i> 27799	Urine	Yellow	Yellow	-	Non-Listeria
<i>Proteus mirabilis</i> 25933	Human vagina	Yellow	Yellow	-	Non-Listeria
<i>Micrococcus luteus</i> 4698	ATCC-undefined	Yellow	Yellow	-	Non-Listeria
<i>Pseudomonas aeruginosa</i> 27853	Blood culture	Yellow	Yellow	-	Non-Listeria
<i>Staphylococcus epidermidis</i> 1228	R-tech labs-undefined	Yellow	Yellow	-	Non-Listeria
<i>Proteus vulgaris</i> 8427	Inner ear infection	Yellow	Yellow	-	Non-Listeria
<i>Streptococcus pyogenes</i> 19615	Pharynx of a child	Yellow	Yellow	-	Non-Listeria

3

4

5

6

7

1 **14.2. Method Comparison Studies**

2

3 Method comparison studies were done for 4 different *Listeria* species paired with four

4 different common food environmental surfaces; *Listeria innocua* (Li) on ceramic tiles,

5 *L.ivanovii* (Liv) on stainless steel, *L. welshmeri* (Lw) on plastic surface, and *L.*


6 *monocytogenes* (Lm) in presence of one log excess *E.coli* on sealed concrete.

7

8

1  
2

Figure 1. Flow chart of method comparison studies



3  
4

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33

#### 14.2.1. Li on Ceramic Tiles:

Serial logarithmic dilutions of overnight culture of Li in TSB were made into sterile TSB. Aliquots containing 100  $\mu$ L of log -6 and log -7 dilutions were plated onto MOX and TSAIE plates for estimating the number of cells loaded onto 4"x 4" ceramic tiles.

Ceramic tiles at 4x4 inches in dimensions were bought from a local Home Depot. They were wrapped in aluminum foil and autoclaved for 15 min at 121 °C. All tiles were kept wrapped at room temperature (RT) until the time of inoculation. The inoculation of surfaces were done at 0.250 mL / surface from the log -4 (high) and log-5 (medium) dilutions, corresponding approximately to 4.67 log CFU and 3.67 log CFU respectively, twenty replicates for each level, per method. Inoculations on surfaces were spread within the corresponding 4"x4" zone with the help of sterile spreaders. Surfaces were left to dry for at least 18 hours at RT.

After the samples were dried on the surfaces for a minimum of 18 hours, they were removed by swabbing using the Tecra Enviroswabs. Both the test method (PDX-LIB) and the reference method (USDA) (10) was done in twenty replicates per level (high and low) of inoculations.

In order to confirm the presence-absence of *Listeria* in all test samples, all tubes of PDX-LIB were streaked onto MOX plates at the end of 48<sup>th</sup> hour of incubation. Dark gray colonies with black zones grown on MOX were forwarded for microscopy and API confirmation. The raw data for Li on tile is given in Appendix I-A. Summary of the data for Li on tile is given in Table 2.

As seen from the Table 2 below, PDX-LIB Labor Saver was more sensitive than the reference method. The reference method resulted in fewer positives than the test method at both levels. Chi squared value also showed that two methods were significantly different from each other at the "medium level" testing of "Li on ceramic tiles".

#### 14.2.2. Liv on Stainless Steel

Serial logarithmic dilutions of overnight culture of Liv in TSB were made into sterile TSB. A one-hundred microliter volume of log -6 and log -7 dilutions was plated onto MOX and TSAIE plates for estimating the number of cells loaded onto 4"x 4" zones on stainless steel.

Stainless steel sheets were bought from Metal Supermarket Inc, Roseville, MN. They were wrapped in aluminum foil and autoclaved for 15 min at 120 °C. Stainless steel sheets were kept wrapped until the time of inoculation. The inoculation of surfaces were done at 0.250 mL / surface from the log -1 (high) and log-2 (medium) dilutions, corresponding approximately to 7.07 log CFU and 6.07 log CFU consecutively, twenty replicates for each level per method. Inoculations on surfaces were spread within the corresponding 4"x4" zone with the help of a disposable sterile spreader. The surfaces were left to dry for a minimum of 18 hours at room temperature.

After the samples were dried on the surfaces for a minimum of 18 hours, they were removed by swabbing using Tecra Enviroswabs. Both test method (PDX-LIB-Labor Saver) and reference method (USDA) (10) was done in twenty replicates per level (high and low) of inoculations.

In order to confirm the presence-absence of *Listeria* in all test samples, all tubes of PDX-LIB Labor Saver were streaked onto MOX plates at the end of 48<sup>th</sup> hour of incubation. Dark gray colonies with black zones grown on MOX were forwarded for microscopy and API confirmation. The raw data for Liv on stainless steel is given in Appendix I-B. The summary of the data for Liv on stainless steel is given in Table 2

As seen from Table 2, PDX-LIB Labor Saver was more sensitive than the reference method. The reference method resulted in far less positive results than the test method at both of the levels tested. Furthermore, chi squared values also suggested that the test method and the reference method were significantly different at both of the levels tested.

#### 14.2.3. Lw on Plastic:

Serial logarithmic dilutions of overnight culture of Lw in TSB were made into sterile TSB. A one-hundred microliter volume of log -6 and log -7 dilutions was plated onto MOX and TSAIE plates for estimating the number of cells loaded onto 4"x 4" zones on plastic.

Polypropylene plastic cutting boards were bought from a local department store. 4"x4" zones were marked on the surface of the cutting boards. They were wrapped in aluminum foil and autoclaved for 15 min at 120 °C. Plastic cutting

1 boards were kept wrapped until the time of inoculation. The inoculation of  
2 surfaces were done at 0.250 mL / surface from the log -4 (high) and log-6 (low)  
3 dilutions, corresponding approximately to 5.18 log CFU and 3.18 log CFU of Lw  
4 consecutively, twenty replicates for each level per method. Inoculations on  
5 surfaces were spread within the corresponding 4"x4" zone with the help of a  
6 disposable sterile spreader. After that surfaces were left to dry min 18 hours at  
7 room temperature.

8  
9 After the samples were dried on the surfaces for a minimum of 18 hours, they  
10 were removed by swabbing using Tecra Enviroswabs. Both test method (PDX-  
11 LIB-Labor Saver) and reference method (USDA) (10) was done in twenty  
12 replicates per level (high and low) of inoculations.

13  
14 In order to confirm the presence of *Listeria* in all test samples, all samples were  
15 streaked onto MOX after at the end of 48<sup>th</sup> hour at 37 °C. Dark gray colonies with  
16 black zones grown on MOX were forwarded for microscopy and API  
17 confirmation. The raw data for Lw on plastic is given in Appendix I-C. The  
18 summary of the data for Lw on plastic is given in Table 2.

19  
20 As seen from Table 2, PDX-LIB Labor Saver was about as sensitive as the  
21 reference method. The reference method gave two more positives at the low level  
22 than the test method, whereas the test method gave one more true positive result  
23 than the reference method at the high level. This difference might be due to the  
24 limited ability of the sampling device to pick up the cells from the rough surface  
25 of the polypropylene cutting board. Chi squared values also suggested that the test  
26 method and the reference method were not significantly different at the low level.

#### 27 28 29 **14.2.4. Lm in 10x *E.coli* on sealed concrete:**

30  
31 Lm in 10x *E.coli* stock culture was made by adding 0.1 mL of overnight grown  
32 Lm ( $10^9$  CFU/mL) and 1 mL of overnight grown *E.coli* ( $10^9$  CFU/mL) into 9  
33 mL of TSB. Serial logarithmic dilutions of Lm in 10x *E.coli* were made into  
34 sterile TSB. A one-hundred microliter volume of log -5, log -6 and log -7 (for *E.*  
35 *coli*) dilutions were plated onto MOX and TSAIE plates for estimating the  
36 number of cells loaded onto 4"x 4" zones on sealed concrete blocks. Number of  
37 cells on MOX was used to estimate number of CFU Lm /mL, and number of  
38 beige colonies on TSAIE was used to estimate the number of CFU *E.coli* / mL.

39  
40 Concrete blocks (15.5"x7.5"x3.5"; length x width x depth) were purchased from a  
41 local Home Depot. Blocks were sealed with a solvent-based concrete sealer (BW  
42 Crete Seal 25 LV, St. Paul, MN). After a minimum of two coats of sealant was  
43 placed on the blocks, they were dried in a chemical hood until all the solvent had  
44 evaporated (minimum of 24 hours). Zones of 4"x4" were marked on the sealed  
45 side using a permanent marker. Before inoculations, sealed concrete blocks were  
46 sprayed with ethanol and allowed to air dry. The ethanol was air dried for at least

1 20 minutes but no more than an hour. The concrete surfaces were not rinsed after  
2 the second application of ethanol was allowed to air dry. The inoculation of sealed  
3 concrete surfaces were done at 0.250 mL / surface from the log -4 *E.coli* / log-5  
4 Lm (high) and log-5 *E.coli* / log-6 Lm (low) dilutions, corresponding  
5 approximately to 3.85 log CFU and 2.85 log CFU respectively, twenty replicates  
6 for each level, per method. Inoculations on surfaces were spread within the  
7 corresponding 4"x4" zone. Surfaces were left to dry for at least 18 hours at RT.  
8

9 After the samples were dried on the surfaces for a minimum of 18 hours, they  
10 were removed by swabbing using Tecra Enviroswabs. Both test method (PDX-  
11 LIB-Labor Saver) and reference method (USDA) (10) was done in twenty  
12 replicates per level (high and low) of inoculations.  
13

14 As seen from Table 2, although the test method yielded in two false negatives in  
15 each level, it performed much better than the reference method. The reference  
16 method missed a total of 12 true positives at the high level and 9 true positives at  
17 the low level. Chi squared values also suggested that there was a significant  
18 difference between the test method and the reference method at both of the levels  
19 tested.  
20

1  
2  
3  
4

Table 2. Summary of the method comparison studies based on 48 hour results

	Levels	Inoculation (CFU/surface)	Replicates	TEST (PDX -LIB Labor Savor)			REFERENCE (USDA)	Chi Squared Value
				Presumptive Positives	Confirmed	False Negatives	Confirmed Positives	
<i>Listeria innocua</i> on ceramic tile	- Ctrl	0	5	0	0	0	0	NA
	Medium	4669	20	19	19	0	6	15.36
	High	46688	20	20	20	0	19	0
<i>Listeria ivanovii</i> on stainless steel	- Ctrl	0	5	0	0	0	0	NA
	Medium	1175000	20	20	20	0	5	20.91
	High	11750000	20	19	19	0	7	13.3
<i>Listeria welshmeri</i> on plastic	- Ctrl	0	5	0	0	0	0	NA
	Low	1506	20	15	15	0	17	1.41
	High	150563	20	20	20	0	19	0
<i>Listeria monocytogenes</i> in 10x <i>E.coli</i> on sealed concrete	- Ctrl	0	5	0	0	0	0	NA
	Medium	715	20	8	9	1	0	7.66
	High	7156	20	14	15	1	3	10.23

1  
2  
3  
4  
5  
6  
7  
8

### **14.3. Ruggedness Studies**

Ruggedness parameters studied were incubation times (28 hours, 32 hours and 46 and 50 hours) and incubation temperatures (36 and 38 °C). Two positive controls (one *L.monocytogenes* and one *L.innocua*) and one negative control (*E.coli*) were tested in 5 replicates at the 1.0E+3 CFU/mL for *Listeria spp* and at 1.0E+7 CFU/mL for *E.coli*. These tests were done on different days as recommended. The summary of the results is given in Table 3.

1  
2  
3  
4

**Table 3. Ruggedness Studies for PDX-LIB**

	Replicates	Variable Incubation hours at fixed incubation temperature (37 °C)						Vairable incubation temperatures at fixed incubation times (30 hours and 48 hours)	
		28 hour	30 hour	32 hour	46hour	48hour	50 hour	36°C*	38°C*
<i>Listeria innocua</i>	A	+	+	+	+	+	+	++	++
	B	+	+	+	+	+	+	++	++
	C	+	+	+	+	+	+	++	++
	D	+	+	+	+	+	+	++	++
	E	+	+	+	+	+	+	++	++
<i>Listeria monocytogenes</i>	A	+	+	+	+	+	+	++	++
	B	+	+	+	+	+	+	++	++
	C	+	+	+	+	+	+	++	++
	D	+	+	+	+	+	+	++	++
	E	+	+	+	+	+	+	++	++
<i>E.coli</i>	A	-	-	-	-	-	-	--	--
	B	-	-	-	-	-	-	--	--
	C	-	-	-	-	-	-	--	--
	D	-	-	-	-	-	-	--	--
	E	-	-	-	-	-	-	--	--
	Neg	-	-	-	-	-	-	-	-

5  
6  
7  
8  
9  
10  
11

**\* Two signs in the each cell represents 30 hour and 48 hour color change profiles;  
+: Black , - : Yellow**

**14.4. Lot to lot variability and shelf stability**

Lot to lot variability studies were done in 5 replicates per level of each microorganism for each of the three production lots tested (Table 4). Shelf life studies were done using three indicator microorganisms in triplicates (Table 5). *L. monocytogenes* ATCC 19119, *L. grayi* (*Listeria grayi* 700545 Lab) and *E. faecalis* ATCC 19433. Five hundred microliters of culture dilutions were inoculated onto the tip of Tecra EnviroSwab, which was then fully submerged into 20 mL of PDX-LIB in its original sterile container.

Table 4: Lot to lot variability studies

	Sample	Lot 50175		Lot 50191		Lot 50200	
		30hr	48hr	30hr	48hr	30hr	48hr
<i>Listeria innocua</i>	Li A	+	+	+	+	+	+
	B	+	+	+	+	+	+
	C	+	+	+	+	+	+
	D	+	+	+	+	+	+
	E	+	+	+	+	+	+
<i>Listeria monocytogenes</i>	Lm A	+	+	+	+	+	+
	B	+	+	+	+	+	+
	C	+	+	+	+	+	+
	D	+	+	+	+	+	+
	E	+	+	+	+	+	+
<i>E. coli</i>	A	-	-	-	-	-	-
	B	-	-	-	-	-	-
	C	-	-	-	-	-	-
	D	-	-	-	-	-	-
	E	-	-	-	-	-	-
	Neg	-	-	-	-	-	-

+: Black , - : Yellow

**Table 5. Shelf life studies of PDX-LIB**

	~ CFU	Lot 63 at 9 months of refrigerated storage (new QC)		
		A	B	C
<i>Listeria monocytogenes</i>	10	+	+	+
	100	+	+	+
<i>Listeria grayi</i>	190	+	+	+
	1900	+	+	+
<i>Enterococcus faecalis</i>	38000	-	-	-
	380000	+	-	-

+: Black , - : Yellow

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46

## 15 Discussion

PDX-LIB Labor Saver is an easy to use and interpret screening test for *Listeria spp* in environmental samples. Inclusivity and exclusivity studies revealed that PDX-LIB Labor Saver is very comprehensive for the detection of *Listeria* species at very low levels (1-52 cells). In regard to specificity, high levels of *Enterococcus* species remain to be a top group for potential false positive results. False positive results, although not desired by the typical end user, still tells a great deal about the overall microbial cleanliness of the areas sampled. Learning about the presence of *Enterococcus* species is important for another aspect, some *Enterococcus* species, namely *E. hirae* produces bacteriocin inhibitory to *Listeria* species (11). Other methods might miss the true presence of *Listeria* if the sample contains *Enterococcus* species.

PDX-LIB Labor Saver was found to be at least as sensitive as the reference method. In fact PDX-LIB Labor Saver has outperformed the reference method in most of the method comparison studies.

Regarding the ruggedness studies, recommended parameters have been studied for PDX-LIB. Results of the ruggedness studies suggested selected deviations from test parameters did not interfere with the true detection of microorganisms selected.

Lot-to-lot variability studies showed that there was no difference between the production lots. Shelf life studies revealed that PDX-LIB is shelf-stable at 10 months of refrigerated storage.

PDX-LIB Labor Saver is a unique, easy to perform rapid detection test for *Listeria* species in environmental samples. Since PDX-LIB Labor Saver is a self-contained test, it minimizes cross contamination. Furthermore PDX-LIB labor saver is compatible with PCR based confirmations (data not shown). This compatibility provides a more economical solution to the *Listeria* screening without any compromise in sensitivity. It is expected that PDX-LIB Labor Saver will revolutionize *Listeria* monitoring in the environmental samples.

## Conclusion

Significant improvement has been made in the PDX-LIB method by simplifying the test protocol. The performance of PDX-LIB Labor Saver has been tested rigorously in AOAC-guideline studies. The results suggested that PDX-LIB Labor Saver was associated with desirable levels of performance characteristics when compared to USDA methods while providing presumptive results in a single stage test.

## 16 References

1. Siragusa, G.R. & Johnson, M. (1990) Monoclonal Antibodies Specific for *Listeria monocytogenes*, *Listeria innocua*, and *Listeria welshimeri*. *Applied and Environmental Microbiology*. **56**:1897-1904.
2. Johnson, J., Jinneman, K., Stelma, G., Smith, B. G., Lye, D., Messer, J., Ulaszek, J., Evsen, L., Gendel, S., Bennett, R. W., Swaminathan, B., Pruckler, J., Steigerwalt, A., Kathariou, S., Yildirim, S., Volokhov, D., Rasooly, a., Chizhikov, V., Wiedmann, M., Fortes, E., Duvall, R. E. & Hitchins, A. D. (2004) Natural Atypical *Listeria innocua* Strains with *Listeria monocytogenes* Pathogenicity Island 1 Gene. *Applied and Environmental Microbiology* **70**:4256-4266.
3. FDA/CFSAN Bad Bug Book, online version at <http://vm.cfsan.fda.gov/~mow/chap6.html>
4. Donnelly, C.W. & Busch, S. (1992) Development of a Repair-Enrichment Broth for Resuscitation of Heat-Injured *Listeria monocytogenes* and *Listeria innocua*. *Applied and Environmental Microbiology*. **58**:14-20.
5. CDC Division of Bacterial and Mycotic Diseases on line publication [http://www.cdc.gov/ncidod/dbmd/diseaseinfo/listeriosis\\_g.htm](http://www.cdc.gov/ncidod/dbmd/diseaseinfo/listeriosis_g.htm)
6. Donnelly, C.W. (2001) *Listeria monocytogenes*: A continuing challenge. *Nutr. Rev.* **59**:183-194.
7. Nogca, H.K., Rudi, K., Naterstad, K., Holck, A.& Lillehaug, D. (2000) Application of 5'-Nuclease PCR for Quantitative Detection of *Listeria monocytogenes* in Pure Cultures, Water, Skim Milk, and Unpasteurized Whole Milk. *Applied and Environmental Microbiology*. **66**:4266-4271.
8. Chen, Y., Ross, W.H., Scott, V.N. & Gombas, D.E. (2003) *Listeria monocytogenes*. Low levels equal low risk. *J. Food Prot.* **66**:570-57.
9. FSIS- Recall Information Center, online publication [http://www.fsis.gov/OA/recalls/rec\\_pr.htm](http://www.fsis.gov/OA/recalls/rec_pr.htm)
10. USDA/FSIS (2002) Microbiology Laboratory Guidelines, Chapter 8; Revision 3; Isolation and Identification of *Listeria monocytogenes* from red meat, poultry, egg, and environmental samples. Online publication in pdf format at <http://www.fsis.usda.gov/OPHS/microlab/mlg8.04.pdf>
11. Siragusa, G.R. (1992) Production of Bacteriocin Inhibitory to *Listeria* species by *Enterococcus hirae*. *Applied and Environmental Microbiology* Volume 58, No. 11, pages 3508-3513.

1

2 **17 Attachments**3 **17.1 Package Inserts: Users Guide and Instruction for Use (IFU) Sheets**

4



5

6

7 **Paradigm Diagnostics *Listeria* Indicator Broth (PDX-LIB)**

8

**User Guide**

9

***Listeria* species Test Kit**

10

**Catalog Number:**

11

**25003-25**

12

**25004-100**

13

**25005-100**

14

**25009-50**

15

16

**Intended Use**

17

18 Paradigm Diagnostics *Listeria* indicator broth (**PDX-LIB**) is intended to be used in the  
19 food processing environment and on food contact surfaces to detect the presence of  
20 *Listeria* species. A color change from yellow to light brown/black is considered  
21 presumptive positive. Applicability of PDX-LIB is limited for selected common *Listeria*  
22 *spp* (*Listeria monocytogenes*, *L. innocua*, *L. ivanovii*, and *L. welshmeri*) on selected  
23 common surface types (Sealed concrete, ceramic tile, stainless steel, and plastic). **AOAC-**  
24 **RI PTM validation studies were conducted at 4 inch x 4-inch surface areas.**

25

26

**Scientific Principle of the Test**

27

28 **PDX-LIB** contains a patented formula of antibiotics, growth enhancers and color  
29 changing compounds. The antibiotics function synergistically to inhibit most non-  
30 *Listeria* microorganisms. Growth enhancers provide recovery nutrients to support the  
31 growth of sub-lethally injured *Listeria*. Indicator compounds will turn the broth from  
32 yellow to black by utilizing the  $\beta$ -glucosidase enzyme produced by *Listeria* species. A  
33 brown to black color 30-48 hours at 37°C indicates a presumptive positive test for  
34 *Listeria spp*. **Positive results can be read as early as 30 hours. Results cannot be**  
35 **considered negative until samples have been incubated for 48 hours.**

36

37

## Additional Notes

Paradigm Diagnostics recommends the use of TECRA EnviroSwab (BioTrace International Bioproducts, Bothell, Washington, 425-398-7993) as sampling device for increased reliability of **PDX-LIB**.

## Materials and Equipment Required

Tecra Enviroswabs for sampling and an incubator capable of maintaining  $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$  (e.g. heat block, water bath, air incubator).

## Confirmation Step

Presumptive positive samples from food contact surfaces or the food processing environment can be confirmed streaking presumptive positive samples onto commonly used selective *Listeria* agar plates such as Modified Oxford Agar or Palcam agars. Typical *Listeria* colonies (dark gray colonies with black zones, generally with dimples) on MOX are used in confirmation protocols given in FDA/BAM (1) or USDA (2) methods.

## Disposal

Decontaminate the **PDX-LIB** after use by autoclave, bleach or other disinfectants in accordance with local, state and federal regulations.

## Product Shelf Life

**PDX-LIB** is stable **for 9 months** at refrigeration temperatures. The expiration date appears on the label along with the lot number. Keep PDX-LIB away from light during storage.

## Precautions

1. *Listeria monocytogenes* is a dangerous human pathogen. When handling samples that possibly contain *L.monocytogenes*, extreme care should be taken to contain the samples and the enriched samples (presumptive positive tubes). Immuno-compromised individuals and pregnant women are particularly endangered by exposure to *L.monocytogenes* and should not be allowed in the vicinity of the testing.
2. WARNING: Some *Enterococci*, particularly *E.hirea* , *E. avium*, *E.feacalis*, and *E. gallinarum* could be resulting in presumptive false positive results. However, all media that noticeably darkens in color must be streaked onto MOX. Because, *Enterococci* do not grow on MOX plates, false positive result remains at a presumptive level.
3. Use of PDX-LIB with some yellow-pigmented sponges should be avoided because of the likely enzyme activity suppression and the resulting possibility of false negative test results.
4. Do not use PDX-LIB past the expiration date that appears on the label
5. Follow standard Good Microbiological Practices where appropriate.

(1) <http://www.cfsan.fda.gov/~ebam/bam-10.html>

(2)[http://www.fsis.usda.gov/Ophs/Microlab/Mlg\\_8\\_04.pdf](http://www.fsis.usda.gov/Ophs/Microlab/Mlg_8_04.pdf)

1 **Warranties and Liabilities**

2  
3 Paradigm Diagnostics Inc warrants the Products manufactured by it will be free from defects in materials  
4 and workmanship when used in accordance with the applicable instructions until the expiration date noted  
5 on the product packaging. Application protocols suggested by Paradigm are intended to be guidelines to  
6 the Buyers of the Products. Each Buyer is expected to validate the applicability of each application  
7 protocol to their individual applications. **PARADIGM DIAGNOSTICS MAKES NO OTHER**  
8 **WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF**  
9 **MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

10  
11 Paradigm Diagnostics' sole obligation with respect to the foregoing warranties shall be, at its option, to  
12 either replace or to refund the purchase price of the Product(s) or part thereof that proves defective in  
13 materials or workmanship within the warranty period, provided the customer notified Paradigm Diagnostics  
14 promptly of any such defect. **PARADIGM DIAGNOSTICS SHALL NOT BE LIABLE FOR ANY**  
15 **DIRECT, INDIRECT OR CONSEQUENTIAL DAMAGES RESULTING FROM ECONOMIC**  
16 **LOSS OR PROPERTY DAMAGES SUSTAINED BY BUYER OR ANY CUSTOMER FROM THE**  
17 **USE OF THE PRODUCT(S).**

18  
19 **For Technical Support, contact Paradigm Diagnostics LLC.**

20  
21  
22 **Paradigm Diagnostics, LLC**  
23 1360 University Ave. W, Suite 455  
24 St. Paul, MN 55104  
25 Web: [www.pdx-inc.com](http://www.pdx-inc.com)  
26 Fax: (612) 545-4657

27  
28 **Technical Support**  
29 Dr. Cem Yurttas  
30 [TechSupport@pdx-inc.com](mailto:TechSupport@pdx-inc.com)  
31 Mobile: 651-226-0381  
32

**Sales Support**  
Frank S. Klisanich  
[SalesSupport@pdx-inc.com](mailto:SalesSupport@pdx-inc.com)  
Mobile: 612-281-0708

33

1  
2 **Instructions for Use: PDX-LIB Labor Saver for Environmental Samples**  
3  
4

5 **1. SAMPLING STEP:** Take an environmental sample following USDA guidelines.

- 6 • Remove recommended sponge out of its sterile wrapping.
- 7 • Sponge the 4"x 4" area to be tested.
- 8 • Return the sponge into the original sampling container.
- 9 • Aseptically add one unit (20 mL) of PDX-LIB on top of sponge in  
10 sampling container, and snap close the sampling container.
- 11 • Incubate upright (sponge fully submerged in PDX-LIB) at 37°C  
12 for 30 to 48 hours.

13  
14 **2. INTERPRETATION STEP:** If color of the media changes from yellow/amber  
15 to brown/black, after **30 to 48** hours of incubation at 37°C, the sample is  
16 considered presumptive positive for *Listeria spp.* **Positive results can be read as**  
17 **early as 30 hours. Results cannot be considered negative until samples have**  
18 **been incubated for 48 hours.** (Note: Although PDX-LIB has been designed to  
19 provide results within 30 hours, 48 hours of incubation at 37 °C has been  
20 shown to significantly improve the sensitivity of the test, and therefore highly  
21 recommended for maximum sensitivity). As with all experimental protocols,  
22 use of negative control (an unused sampling device containing one unit of PDX-  
23 LIB, incubated along side with the environmental samples) in each set of samples  
24 is recommended.

25  
26 **NOTE: For audio-visual instruction, please visit us at**

27 <http://www.pdx-inc.com/training/index.html>



28 **and click on the link below**

29 [Training Video - PDX-LIB-List 'Labor Saver' Method](#)  
30  
31  
32

- 1
- 2
- 3
- 4



THE EASIEST SOLUTION TO ENVIRONMENTAL *LISTERIA* MONITORING  
PDX-LIB LABOR SAVER



For more Information: Call Paradigm Diagnostics at 612-281-0708

1  
2  
3  
4

**APPENDIXES**

**Appendix I. Method Comparison Data A) Li on Tile**

Medium	PDX-LIB				USDA			
	30 hour	48 hour	MOX	API	24 hour	FB	FB MOX	API
A	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
B	-	+	+	<i>Listeria innocua</i>	+	+		
C	-	-	-		-	+	+	<i>Listeria innocua</i>
D	+	+	+	<i>Listeria innocua</i>	+	+		
E	+	+	+	<i>Listeria innocua</i>	-	-		
F	-	+	+	<i>Listeria innocua</i>	-	-		
G	-	+	+	<i>Listeria innocua</i>	-	-		
H	-	+	+	<i>Listeria innocua</i>	-	-		
I	+	+	+	<i>Listeria innocua</i>	-	-		
J	+	+	+	<i>Listeria innocua</i>	-	-		
K	-	+	+	<i>Listeria innocua</i>	-	+	+	<i>Listeria innocua</i>
L	+	+	+	<i>Listeria innocua</i>	-	-		
M	+	+	+	<i>Listeria innocua</i>	-	+	+	<i>Listeria innocua</i>
N	+	+	+	<i>Listeria innocua</i>	-	-		
O	+	+	+	<i>Listeria innocua</i>	-	-		
P	+	+	+	<i>Listeria innocua</i>	-	+	+	<i>Listeria innocua</i>
Q	+	+	+	<i>Listeria innocua</i>	-	-		
R	-	+	+	<i>Listeria innocua</i>	-	-		
S	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
T	-	+	+	<i>Listeria innocua</i>	-	-		
High Sample	PDX-LIB				USDA			
Sample	30 hour	48 hour	MOX	API	24 hour	FB	FB MOX	API
A	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
B	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
C	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
D	+	+	+	<i>Listeria innocua</i>	-	+	+	<i>Listeria innocua</i>
E	+	+	+	<i>Listeria innocua</i>	-	+	+	<i>Listeria innocua</i>
F	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
G	+	+	+	<i>Listeria innocua</i>	-	-		
H	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
I	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
J	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
K	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
L	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
M	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
N	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
O	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
P	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
Q	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
R	+	+	+	<i>Listeria innocua</i>	+	+	+	<i>Listeria innocua</i>
S	+	+	+	<i>Listeria innocua</i>	-	+		<i>Listeria innocua</i>
T	+	+	+	<i>Listeria innocua</i>	+	+		<i>Listeria innocua</i>
Blk 1	-	-	-		-	-		
Blk 2	-	-	-		-	-		
Blk 3	-	-	-		-	-		
Blk 4	-	-	-		-	-		
Blk 5	-	-	-		-	-		

5

6  
7  
8  
9  
10

24 hour=Presence of *Listeria* on MOX from 24 hour UVM; FB: Darkening of Fraser Broth; FB MOX: presence of *Listeria* colonies on MOX from darkened FB.

1

**Appendix I. Method Comparison Data B) Liv on stainless steel**

Medium	PDX-LIB				USDA			
	30 hour	48 hour	MOX	API	24 hour	FB	FB MOX	API
A	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
B	-	+	+	<i>Listeria ivanovii</i>	-	-		
C	-	+	+	<i>Listeria ivanovii</i>	-	-		
D	-	+	+	<i>Listeria ivanovii</i>	-	-		
E	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
F	-	+	+	<i>Listeria ivanovii</i>	-	-		
G	-	+	+	<i>Listeria ivanovii</i>	-	-		
H	-	+	+	<i>Listeria ivanovii</i>	-	-		
I	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
J	-	+	+	<i>Listeria ivanovii</i>	-	-		
K	-	+	+	<i>Listeria ivanovii</i>	-	-		
L	-	+	+	<i>Listeria ivanovii</i>	-	-		
M	-	+	+	<i>Listeria ivanovii</i>	-	-		
N	-	+	+	<i>Listeria ivanovii</i>	-	-		
O	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
P	-	+	+	<i>Listeria ivanovii</i>	-	-		
Q	-	+	+	<i>Listeria ivanovii</i>	-	-		
R	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
S	-	+	+	<i>Listeria ivanovii</i>	-	-		
T	-	+	+	<i>Listeria ivanovii</i>	-	-		
High	PDX-LIB				USDA			
Sample	30 hour	48 hour	MOX	API	24 hour	FB	FB MOX	API
A	-	+	+	<i>Listeria ivanovii</i>	-	-		
B	-	+	+	<i>Listeria ivanovii</i>	-	-		
C	-	-	-		-	-		
D	-	+	+	<i>Listeria ivanovii</i>	-	-		
E	-	+	+	<i>Listeria ivanovii</i>	-	-		
F	+	+	+	<i>Listeria ivanovii</i>	-	-		
G	-	+	+	<i>Listeria ivanovii</i>	-	-		
H	-	+	+	<i>Listeria ivanovii</i>	-	-		
I	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
J	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
K	+	+	+	<i>Listeria ivanovii</i>	-	-		
L	+	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
M	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
N	-	+	+	<i>Listeria ivanovii</i>	-	-		
O	-	+	+	<i>Listeria ivanovii</i>	-	-		
P	+	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
Q	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
R	+	+	+	<i>Listeria ivanovii</i>	-	-		
S	-	+	+	<i>Listeria ivanovii</i>	-	+	+	<i>Listeria ivanovii</i>
T	-	+	+	<i>Listeria ivanovii</i>	-	-		
Blk 1	-	-	-		-	-		
Blk 2	-	-	-		-	-		
Blk 3	-	-	-		-	-		
Blk 4	-	-	-		-	-		
Blk 5	-	-	-		-	-		

2

3  
4  
5  
6  
7

24 hour=Presence of *Listeria* on MOX from 24 hour UVM; FB: Darkening of Fraser Broth; FB MOX: presence of *Listeria* colonies on MOX from darkened FB.

1  
2  
3

**Appendix I. Method Comparison Data C) Lw on plastic**

Low	PDX-LIB				USDA			
	30 hour	48 hour	MOX	API	24 hour	FB	FB MOX	API
A	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
B	-	+	+	<i>Listeria welshmeri</i>	-	+	+	<i>Listeria welshmeri</i>
C	-	+	+	<i>Listeria welshmeri</i>	-	+	+	<i>Listeria welshmeri</i>
D	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
E	-	+	+	<i>Listeria welshmeri</i>	-	+	+	<i>Listeria welshmeri</i>
F	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
G	-	+	+	<i>Listeria welshmeri</i>	-	-		
H	-	-	-		-	+		<i>Listeria welshmeri</i>
I	-	-	-		+	+		<i>Listeria welshmeri</i>
J	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
K	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
L	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
M	-	+	+	<i>Listeria welshmeri</i>	-	-		
N	-	+	+	<i>Listeria welshmeri</i>	-	-		
O	-	-	-		+	+		<i>Listeria welshmeri</i>
P	-	-	-		+	+		<i>Listeria welshmeri</i>
Q	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
R	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
S	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
T	-	-	-		+	+		<i>Listeria welshmeri</i>
High	PDX-LIB				USDA			
Sample	30 hour	48 hour	MOX	API	24 hour	FB	FB MOX	API
A	+	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
B	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
C	+	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
D	+	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
E	+	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
F	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
G	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
H	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
I	+	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
J	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
K	-	+	+	<i>Listeria welshmeri</i>	-	-		
L	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
M	+	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
N	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
O	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
P	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
Q	+	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
R	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
S	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
T	-	+	+	<i>Listeria welshmeri</i>	+	+		<i>Listeria welshmeri</i>
Blk 1	-	-	-		-	-		
Blk 2	-	-	-		-	-		
Blk 3	-	-	-		-	-		
Blk 4	-	-	-		-	-		
Blk 5	-	-	-		-	-		

4

24 hour=Presence of *Listeria* on MOX from 24 hour UVM; FB: Darkening of Fraser Broth; FB MOX: presence of *Listeria* colonies on MOX from darkened FB.

5  
6  
7  
8  
9  
10  
11

1  
2  
3

**Appendix I. Method Comparison Data D) Lm in 10x *E.coli* on sealed concrete**

Medium	PDX-LIB				USDA			
	30 hour	48 hour	MOX	API	24 hour	FB	FB MOX	API
A	-	-	-		-	-		
B	-	+	+	<i>Listeria monocytogenes</i>	-	-		
C	-	+	+	<i>Listeria monocytogenes</i>	-	-		
D	-	-	-		-	-		
E	-	-	-		-	-		
F	-	-	-		-	-		
G	-	+	+	<i>Listeria monocytogenes</i>	-	-		
H	-	-	-		-	-		
I	-	-	-		-	-		
J	-	+	+	<i>Listeria monocytogenes</i>	-	-		
K	-	+	+	<i>Listeria monocytogenes</i>	-	-		
L	-	+	+	<i>Listeria monocytogenes</i>	-	-		
M	-	-	-		-	-		
N	-	-	-		-	-		
O	-	+	+	<i>Listeria monocytogenes</i>	-	-		
P	-	-	-		-	-		
Q	+	+	+	<i>Listeria monocytogenes</i>	-	-		
R	-	-	-		-	-		
S	-	-	+	<i>Listeria monocytogenes</i>	-	-		
T	-	-	-		-	-		
High	PDX-LIB				USDA			
Sample	30 hour	48 hour	MOX	API	24 hour	FB	FB MOX	API
A	-	-	-		-	-		
B	+	+	+	<i>Listeria monocytogenes</i>	-	-		
C	-	+	+	<i>Listeria monocytogenes</i>	-	-		
D	+	+	+	<i>Listeria monocytogenes</i>	-	-		
E	-	+	+	<i>Listeria monocytogenes</i>	-	-		
F	+	+	+	<i>Listeria monocytogenes</i>	-	-		
G	-	+	+	<i>Listeria monocytogenes</i>	+	-		<i>Listeria monocytogenes</i>
H	-	+	+	<i>Listeria monocytogenes</i>	-	-		
I	-	+	+	<i>Listeria monocytogenes</i>	-	-		
J	-	-	-		-	-		
K	-	+	+	<i>Listeria monocytogenes</i>	-	+	+	<i>Listeria monocytogenes</i>
L	-	+	+	<i>Listeria monocytogenes</i>	-	-		
M	-	+	+	<i>Listeria monocytogenes</i>	-	-		
N	-	+	+	<i>Listeria monocytogenes</i>	+	+		<i>Listeria monocytogenes</i>
O	-	-	-		-	-		
P	-	-	-		-	-		
Q	-	+	+	<i>Listeria monocytogenes</i>	-	-		
R	-	+	+	<i>Listeria monocytogenes</i>	-	-		
S	-	-	-		-	-		
T	-	-	+	<i>Listeria monocytogenes</i>	-	-		
Blk 1	-	-	-		-	-		
Blk 2	-	-	-		-	-		
Blk 3	-	-	-		-	-		
Blk 4	-	-	-		-	-		
Blk 5	-	-	-		-	-		

4

5  
6  
7  
8

24 hour=Presence of *Listeria* on MOX from 24 hour UVM; FB: Darkening of Fraser Broth; FB MOX: presence of *Listeria* colonies on MOX from darkened FB.